

SteriGraft™ – Tibia Whole

Product Code	Description	Size & Measurements
490	Tibia Whole	Condyle to Condyle

Source and Common Surgical Applications: Same as below

Preservation Method: Frozen
Containers must be stored at -40°C or colder



SteriGraft™ – Tibia Section

Product Code	Description	Size & Measurements
406	Tibia Section, 80 – 100 mm	L = 80 - 100 mm D = mm (varies)
407	Tibia Section, 101 - 200 mm	L = 101 - 200 mm D = mm (varies)

Source: Tibia

Preservation Method: Lyophilized
Containers must be stored at ambient temperature

Common Surgical Applications: Spinal fusion or general orthopaedic reconstruction



SteriGraft™ – Tibia Section - Frozen

Product Code	Description	Size & Measurements
494	Tibia Section, 100 - 200 mm - Frozen	L = 100 - 200 mm D = mm (varies) ID = mm (varies)

Source and Common Surgical Applications: Same as above

Preservation Method: Frozen
Containers must be stored at -40°C or colder

