

SteriGraft™ – Femur Section


Product Code	Description	Size & Measurements
604	Femur Section, 100 mm	L = 98 - 102 mm D ≥ 20 mm
605	Femur Section, 150 mm	L = 148 - 152 mm D ≥ 20 mm
606	Femur Section, 200 mm	L = 198 - 202 mm D ≥ 20 mm
609	Femur Section, 250 mm	L = 248 - 252 mm D ≥ 20 mm



Source: Femur

Preservation Method: Lyophilized
Containers must be stored at ambient temperature

Common Surgical Applications: Lumbar fusion, hip revision or general orthopaedic reconstruction



SteriGraft™ – Femur Section (Frozen)

Product Code	Description	Size & Measurements
699	Femur Section, 150 mm - Frozen	L = 148 - 152 mm D ≥ 20 mm
697	Femur Section, 200 mm - Frozen	L = 198 - 202 mm D ≥ 20 mm
694	Femur Section, 250 mm - Frozen	L = 248 - 252 mm D ≥ 20 mm

Source & Common Surgical Applications: Same as above

Preservation Method: Frozen
Containers must be stored at -40°C or colder