


SteriGraftTM – Cortical Ring Castled and Lordosed 5°


Product Code	Description	Size & Measurements
70005	Cortical Ring Castled and Lordosed 5°, 5 mm	H = 5 mm W = 13 mm L = 11 mm
70006	Cortical Ring Castled and Lordosed 5°, 6 mm	H = 6 mm W = 13 mm L = 11 mm
70007	Cortical Ring Castled and Lordosed 5°, 7 mm	H = 7 mm W = 13 mm L = 11 mm
70008	Cortical Ring Castled and Lordosed 5°, 8 mm	H = 8 mm W = 13 mm L = 11 mm
70009	Cortical Ring Castled and Lordosed 5°, 9 mm	H = 9 mm W = 13 mm L = 11 mm
70010	Cortical Ring Castled & Lordosed 5°, 10 mm	H = 10 mm W = 13 mm L = 11 mm
70011	Cortical Ring Castled & Lordosed 5°, 11 mm	H = 11 mm W = 13 mm L = 11 mm



Source: Radius or ulna

Preservation Method: Lyophilized
Containers must be stored at ambient temperature

Common Surgical Applications: Cervical spinal fusion



SteriGraftTM – Cortical Ring Castled

Product Code	Description	Size & Measurements
71005	Cortical Ring Castled, 5 mm	H: 5mm (+0.3/-0.2) W: 13mm (+/-1.0mm) L: 11mm (+/-1.0mm)
71006	Cortical Ring Castled, 6 mm	H: 6mm (+0.3/-0.2) W: 13mm (+/-1.0mm) L: 11mm (+/-1.0mm)
71007	Cortical Ring Castled, 7 mm	H: 7mm (+0.3/-0.2) W: 13mm (+/-1.0mm) L: 11mm (+/-1.0mm)
71008	Cortical Ring Castled, 8 mm	H: 8mm (+0.3/-0.2) W: 13mm (+/-1.0mm) L: 11mm (+/-1.0mm)
71009	Cortical Ring Castled, 9 mm	H: 9mm (+0.3/-0.2) W: 13mm (+/-1.0mm) L: 11mm (+/-1.0mm)
71010	Cortical Ring Castled, 10 mm	H: 10mm (+0.3/-0.2) W: 13mm (+/-1.0mm) L: 11mm (+/-1.0mm)
71011	Cortical Ring Castled, 11 mm	H: 11mm (+0.3/-0.2) W: 13mm (+/-1.0mm) L: 11mm (+/-1.0mm)