


SteriGraftTM – Unicortical Cancellous Lordotic Bone Wedge

Product Code	Description	Size & Measurements
340	Unicortical Cancellous Lordotic Bone Wedge (Lg)	H = 12 mm L = 38 mm W = 22 mm (+/-1.0mm)
341	Unicortical Cancellous Lordotic Bone Wedge (Sm)	H = 6 mm L = 24 mm W = 14 mm (+/-1.0mm)




Source: Ilium, talus, calcaneus, femoral condyles, proximal tibia, or patella

Preservation Method: Lyophilized
Containers must be stored at ambient temperature

Common Surgical Applications: Foot / Ankle Reconstruction

SteriGraftTM – Lordotic Bone Wedge

Product Code	Description	Size & Measurements
330	Lordotic Bone Wedge (Lg)	H = 10 mm L = 50 mm W = 40 mm (+/-0.5 mm)
331	Lordotic Bone Wedge (Sm)	H = 6 mm L = 24 mm W = 14 mm (+/-0.5 mm)



Source: Ilium, talus, calcaneus, femoral condyles, proximal tibia, or patella

Preservation Method: Lyophilized
Containers must be stored at ambient temperature

Common Surgical Applications: Foot / Ankle Reconstruction