


SteriGraft™ – Cortical Plate

Product Code	Description	Size & Measurements
440	Cortical Plate 10 x 100 mm	L = 98 - 102 mm W = 9 - 11 mm T ≥ 3 mm
441	Cortical Plate 10 x 150 mm	L = 148 - 152 mm W = 9 - 11 mm T ≥ 3 mm
442	Cortical Plate 10 x 200 mm	L = 198 - 202 mm W = 9 - 11 mm T ≥ 3 mm
443	Cortical Plate 20 x 100 mm	L = 98 - 102 mm W = 19 - 21 mm T ≥ 3 mm
444	Cortical Plate 20 x 150 mm	L = 148 - 152 mm W = 19 - 21 mm T ≥ 3 mm
445	Cortical Plate 20 x 200 mm	L = 198 - 202 mm W = 19 - 21 mm T ≥ 3 mm



Source: Femur or tibia

Preservation Method: Lyophilized
Containers must be stored at ambient temperature

Common Surgical Applications: Hip revision or general orthopaedic reconstruction

SteriGraft™ – Cortical Plate – Frozen

Product Code	Description	Size & Measurements
599	Cortical Plate 20 x 200 mm - Frozen	L = 198 - 202 mm W = 19 - 21 mm T ≥ 3 mm

Preservation Method: Frozen
Containers must be stored at -40°C or colder

Source & Common Surgical Applications: Same as above