

**SteriGraft<sup>TM</sup> – Humerus - Frozen**


Product Code	Description	Size & Measurements
291	Proximal Humerus Right – Frozen	L ≥ 160 mm D = mm (varies) D2 = mm (varies) See diagram
292	Proximal Humerus Left - Frozen	L ≥ 160 mm D = mm (varies) D2 = mm (varies) See diagram
298	Distal Humerus RT/LT - Frozen	L ≥ 150 mm W' = mm (varies) W = mm (varies) D' = mm (varies)



Source: Humerus

Preservation Method: Frozen  
*Containers must be stored at -40°C or colder*

Common Surgical Applications: ORIF or total shoulder reconstruction




*Drawing for 291 / 292*

**SteriGraft<sup>TM</sup> – Humerus Section – Frozen**

Product Code	Description	Size & Measurements
293	Humerus Section - Frozen	L ≥ 148 mm; D ≥ 15 mm D2 = mm (varies)

Source & Common Surgical Applications: Same as above.

Preservation Method: Frozen  
*Containers must be stored at -40°C or colder*




*Drawings for 298*

**SteriGraft<sup>TM</sup> – Humerus Section – Freeze Dried**

Product Code	Description	Size & Measurements
923	Humerus Section	L ≥ 148 mm; D ≥ 15 mm D2 = mm (varies)

Source & Common Surgical Applications: Same as above.

Preservation Method: Lyophilized  
*Containers must be stored at ambient temperature*



*Drawing for 293 / 923*